

fMRI Lab – Pre-processing and GLM modeling

This is an example of blocked fMRI data acquired on a modified 2T Siemens Magnetom Vision System. There are 96 functional scans (BOLD/EPI images; fM00223_004 – fM00223_009 and one structural T1-weighted image (sM00223_002).

The blocked design task is auditory stimulation (two-syllable words presented at a rate of 1 per second) alternating with periods of silence/rest.

The BOLD/EPI parameters are:


- 64 contiguous slices (64x64x64)
- Voxel size: 3mm x 3mm x 3mm
- TR = 7s


In this lab, you will:

- 1) normalize the structural and functional data
- 2) smooth the functional data
 - You can play around with different smoothing kernels or examine the difference between smoothing and not smoothing the data
- 3) specify a design matrix (model)
- 4) estimate the model
- 5) inspect the results
- 6) specify contrasts.



I. Launch SPM12 from Matlab

II. Co-register structural and functional data:





- A. For the “Reference Image,” select sM00223_002.img and then click “Done”
- B. For the “Source Image,” select the first functional image in the functional image directory: “fM00223_004.img” and then click “Done”



When all the files have been specified in the Batch Editor, click on the green arrow.



When the co-registration is done, you will see this in the SPM Graphics window:

Normalised Mutual Information Coregistration


$X1 = 3.000^{\circ}X + 0.018^{\circ}Y + 0.017^{\circ}Z + 28.792$

$Y1 = -0.018^{\circ}X + 3.000^{\circ}Y + 0.042^{\circ}Z + 32.813$

$Z1 = -0.008^{\circ}X - 0.014^{\circ}Y + 1.000^{\circ}Z - 8.475$




II. Segment the structural data





- 1) Click on Volumes
- 2) Click Specify
- 3) Navigate to the sM00223 directory and highlight "M00223_002.img"; then click "Done"


- 4) After you see the structural image in the Volumes field, click the green button to segment. This will take a few of minutes.
- 5) When the segmentation is complete, in the sM00223 directory, you should see 6 new files: c1sM00223_002, c2sM00223_002, c3sM00223_002, c4_sM00223_002, c5sM00223_002, and sM00223_002_seg8.




III. Normalize the structural data:



- 1) After clicking "Specify", highlight "Image to Align" and select c1sM00223_002.nii (click Done).




2) Highlight “Images to Write”; select “sM00223_002.img”; click “Done”:





3) After selecting both files, click on the green arrow to Normalize. This step will take several minutes.

- a. When this step is finished, there will be 3 new files in the sM00223 folder:
 - wsM00223_002.img,
 - wsM00223_002.hdr, (normalized structural image) and
 - y_cs1sM00223_002 (deformation field)




IV. Apply the deformation field to the co-registered functional data:



- 1) choose "Normalize (Write)"
- 2) click "Data" and "specify":



- 3) For "Deformation Field," select "y_c1sM00223_002.nii" in the structural folder and click "Done"
- 4) For "Images to write," select all the functional images in the "fM00223" folder then click the green arrow.
- 5) After this step is complete, there will be 192 new functional images with the prefix "w" (e.g., wfM00223_004).
- 6) Move these 'normalized' functional images into a separate folder called Normalized within the "fM00223" folder




V. Smooth the normalized functional data.




- 1) For “Images to Smooth,” select all the normalized functional scans and then click done.
- 2) Then click the green arrow.
 - a. You can leave the FWHM at the default 8,8,8 for now. If there is time, you can play around with different smoothing kernels or not smoothing the functional data at all.
- 3) When this step is finished, you should have 192 new files with the prefix “s” in the Normalized folder (e.g., swfM00223_004).
- 4) Move these smoothed, normalized functional images into a separate folder called Smooth888 within the “fM00223” folder
- 5) Because there were no “dummy” lead in scans in the experiment, we can remove the first few scans (i.e., you can move the first 6 scans into a separate folder so the experiment starts with a sound block).

VI. Specify the design matrix:

- 1) Within the MoAEpilot directory, make a directory called “Stats888” to put the files generated from the 1st level analysis.
- 2) Click on “Specify 1st-level”





- 3) For Directory, select the Stats888 folder you just created
- 4) Units for design, select seconds
- 5) Interscan interval, 7
- 6) For Data & Design, select New Subject/Session
- 7) For Scans, select all smoothed, normalized, functional scans
- 8) Click on Conditions and select "New: Condition"
 - a. Name, call first condition something like "Sound"
 - b. Onsets, this is asking for the time in seconds when the sound blocks occur.

VII. Estimated the model:


- 1) Click Estimate
- 2) Click Select SPM.mat, then specify
- 3) Click the SPM.mat file in the Stats888 folder
- 4) Click green arrow.



VII. Estimated the model:

- 1) Click Results
- 2) Select the SPM.mat file in the Stats888 folder

VII. Estimated the model:

- 1) Select t-contrasts
- 2) Then click "Define new contrast"






Contrast 1: Sound > Rest, specify 1 -1


Contrast 2: Rest > sound, specify -1 1

To view results:

- 1) Click results
- 2) Click SPM.mat file, done
- 3) Pick contrast to look at, select 001{T} Sound > rest, then Done
- 4) Apply masking: none
- 5) p value adjustment to control: FEW
- 6) p value (PWE): 0.05
- 7) & extent threshold [voxels]: 0
- 8) You should see Glass Brain rendition of activations



9) You can also over lay the activations on the subject's structural image:



- 10) Can overlay activations on slices or section (this option allows you to click around a 3D volume)
- Select normalize structural image for the overlay.

Auditory Experiment Timing Parameters

fM00223_004	1 Rest (discard)		Onsets
fM00223_005			
fM00223_006			
fM00223_007			
fM00223_008			
fM00223_009			
		seconds	
fM00223_010	2 Sound	1-7	
fM00223_011		8-14	
fM00223_012		15-21	
fM00223_013		22-28	
fM00223_014		29-35	
fM00223_015		36-42	
fM00223_016	3 Rest	43-49	
fM00223_017		50-56	
fM00223_018		57-63	
fM00223_019		64-70	
fM00223_020		71-77	
fM00223_021		78-84	
fM00223_022	4 Sound	85-91	
fM00223_023		92-98	
fM00223_024		99-105	
fM00223_025		106-112	
fM00223_026		113-119	
fM00223_027		120-126	
fM00223_028	5 Rest	127-133	
fM00223_029		134-140	
fM00223_030		141-147	
fM00223_031		148-154	
fM00223_032		155-161	
fM00223_033		162-168	
fM00223_034	6 Sound	169-175	
fM00223_035		176-182	
fM00223_036		183-189	
fM00223_037		190-196	
fM00223_038		197-203	
fM00223_039		204-210	
fM00223_040	7 Rest	211-217	
fM00223_041		218-224	
fM00223_042		225-231	
fM00223_043		232-238	
fM00223_044		239-245	
fM00223_045		246-252	
fM00223_046	8 Sound	253-259	
fM00223_047		260-266	
fM00223_048		267-273	
fM00223_049		274-280	
fM00223_050		281-287	

fM00223_051		288-294
fM00223_052	9	295-301
fM00223_053	Rest	302-308
fM00223_054		309-315
fM00223_055		316-322
fM00223_056		323-329
fM00223_057		330-336
fM00223_058	10	337-343
fM00223_059	Sound	344-350
fM00223_060		351-357
fM00223_061		358-364
fM00223_062		365-371
fM00223_063		372-378
fM00223_064	11	379-385
fM00223_065	Rest	386-392
fM00223_066		393-399
fM00223_067		400-406
fM00223_068		407-413
fM00223_069		414-420
fM00223_070	12	421-427
fM00223_071	Sound	428-434
fM00223_072		435-441
fM00223_073		442-448
fM00223_074		449-455
fM00223_075		456-462
fM00223_076	13	463-469
fM00223_077	Rest	470-476
fM00223_078		477-483
fM00223_079		484-490
fM00223_080		491-497
fM00223_081		498-504
fM00223_082	14	505-511
fM00223_083	Sound	512-518
fM00223_084		519-525
fM00223_085		526-532
fM00223_086		533-539
fM00223_087		540-546
fM00223_088	15	547-553
fM00223_089	Rest	554-560
fM00223_090		561-567
fM00223_091		568-574
fM00223_092		575-581
fM00223_093		582-588
fM00223_094	16	589-595
fM00223_095	Sound	596-602
fM00223_096		603-609
fM00223_097		610-616
fM00223_098		617-623
fM00223_099		624-630